FUNCTIONAL AND STRUCTURAL IMAGING WITH MICROBUBBLE CONTRAST AGENTS



robert.eckersley@kcl.ac.uk

Imaging Sciences Division

Biomedical Engineering Department King's College London England, U.K.



Overview

- Methods to improve quantitative measure of microbubble contrast:
 - attenuation correction
 - non-linear artifact correction
- A different approach to bubble counting.
- Method to reveal fine structural detail in the vasculature (beyond the diffraction limited system resolution).

Tang, M. X., et al. (2011). Quantitative contrastenhanced ultrasound imaging: a review of sources of variability. *Interface Focus*, 1(4), 520-539.

Attenuation Compensation



Assume uniform intensity in lumen.


Remove effects of speckle

Rescale and normalise tissue region below lumen

Attenuation correction and normalisation for quantification of contrast enhancement in ultrasound images of carotid arteries.

Cheung et al


(ACCEPTED in Ultrasound in Medicine and Biology 21/02/2015)



Nonlinear propagation artefact in contrast US



Courtesy of Prof. David Cosgrove



Courtesy of Prof. Roxy Senior


 Nonlinear propagation causes qualitative (misclassification of tissue) and quantitative (error in concentration estimation) artefacts

Tang, M.-X., N. Kamiyama, and R.J. Eckersley, Effects of Nonlinear Propagation in Ultrasound Contrast Agent Imaging.

Ultrasound in Medicine & Biology, 2010. **36**(3): p. 459-466

Tang, M.X. and R.J. Eckersley, Nonlinear propagation of ultrasound through microbubble contrast agents and implications for Imaging. IEEE Transactions on Ultrasonics Ferroelectrics and Frequency Control, 2006. **53**(12): p. 2406-2415

Correction of non-linear propagation: *In-vitro* demonstration



In-Vivo Velocity Mapping




Fig 1. Super-resolution image of mouse ear microvasculature


Fig 2. Image of the directional component of velocity.


Fig 3. Image of the speed of bubble flow.

Correction of non-linear propagation: In-vivo demonstration

Estimated tissue artefact image











Correction of nonlinear propagation artifact in contrast enhanced ultrasound imaging of carotid arteries: methods and in vitro evaluation. Yildiz et al 2015.

(ACCEPTED 11th March 2015 for publication in Ultrasound in Medicine and Biology)













TSD






TSD









Corrected PI image

In vitro results







Demodulated RF data averaged over ten central columns in the original versus corrected PI image Contrast to Tissue Ratio(CTR) for the original and corrected images. CTR increases significantly after correction.

Correction of non-linear propagation: In-vivo demonstration

Estimated tissue artefact image



Corrected PI image



Correction of nonlinear propagation artifact in contrast enhanced ultrasound imaging of carotid arteries: methods and in vitro evaluation. Yildiz et al 2015.

(ACCEPTED 11th March 2015 for publication in Ultrasound in Medicine and Biology)