

Ultrasound in Med. & Biol., Vol. 40, No. 1, pp. 168–176, 2014
Copyright © 2014 World Federation for Ultrasound in Medicine & Biology
Printed in the USA. All rights reserved
0301-5629/\$ - see front matter

<http://dx.doi.org/10.1016/j.ultrasmedbio.2013.08.020>

● *Original Contribution*

DEVELOPMENT OF OIL-IN-GELATIN PHANTOMS FOR VISCOELASTICITY MEASUREMENT IN ULTRASOUND SHEAR WAVE ELASTOGRAPHY

MAN M. NGUYEN,^{*†} SHIWEI ZHOU,^{*} JEAN-LUC ROBERT,^{*} VIJAY SHAMDASANI,[‡] and HUA XIE^{*}

^{*}Philips Research North America, Briarcliff Manor, New York, USA; [†]University of Southern California, University Park Campus, Los Angeles, California, USA; and [‡]Philips Healthcare, Bothell, Washington, USA

Hua Xie
Philips Research North America
2014-02-21

PHILIPS


Study Motivations and Objectives

Work was carried out in 2012 with Dr. Man Nguyen (summer intern)

- Develop and fabricate viscous tissue mimicking phantoms for ultrasound shear wave elastography
- Develop and evaluate processing and reconstruction algorithms that construct elasticity and viscosity
- Clinical applications of liver tissue viscosity measurement
 - This study was mainly motivated by Barry et al. “Shear wave dispersion measures liver steatosis”, *Ultrasound Med Biol* 2012;38(2):175-182
 - Some clinical studies reported higher liver stiffness and correlated increased shear viscosity in patients with fibrotic liver diseases

PHILIPS

Methods



Viscous tissue-mimicking phantoms

- 3-7% Gelatin (elasticity)
- 0-40% Castor oil (elasticity & viscosity)
- 2% graphite (backscattering)
- 7% propanol (Speed of Sound)

3

Methods: Material Concentrations of Viscous Phantoms



Table 1. Material concentrations by weight for oil-in-gelatin phantoms.

Material Components	Phantoms				
	1-4	5-6	7-8	9-12	13- 14
Gelatin (%)	3	4	5	6	7
Graphite (%)	2	2	2	2	2
N-propanol (%)	7	7	7	7	7
Castor oil (%)	0, 10, 20, 40	0, 20	0, 20	0, 10, 20, 40	0, 20

Example: material weight for 4% and 7% gelatin phantoms

Phantom	Gelatin	Castor Oil	DI water	(2-propanol)	Graphite	Total weight
4% phantoms	4% (40g)	0% (0g)	87% (870g)	7% (70g)	2% (20g)	100% (1000g)
	4% (40g)	20% (200g)	67% (670g)	7% (70g)	2% (20g)	100% (1000g)
7% phantoms	7% (70g)	0% (0g)	84% (840g)	7% (70g)	2% (20g)	100% (1000g)
	7% (70g)	20% (200g)	64% (640g)	7% (70g)	2% (20g)	100% (1000g)

PHILIPS



Methods: Viscoelasticity Measurements

- Shear wave elastography research prototype on Philips iU22 & C5-1 curvilinear transducer
- Push and tracking pulse specs
 - Push pulses with pulse repetition frequency of 70 Hz → shear waves with fundamental at 70 Hz and harmonics at 140, 210, 280 Hz
 - Tracking PRF 2KHz at multiple tracking locations
- Data processing steps
 - Motion tracking to estimate shear wave induced displacement
 - Bandpass filter to remove background motion
 - Shear wave phase velocity measurement ([70, 280 Hz])
 - Voigt model fitting


$$c_s(\omega_s) = \sqrt{\frac{2(\mu_1^2 + \omega_s^2 \mu_2^2)}{\rho(\mu_1 + \sqrt{\mu_1^2 + \omega_s^2 \mu_2^2})}}$$

c_s : shear wave speed
 μ_1 : shear modulus
 μ_2 : shear viscosity
 ρ : medium density


7

Acoustic Properties of Phantoms as a Function of Oil Concentration

Speed of Sound




Attenuation Coefficient



- Fabricated phantoms have comparable sound speeds and attenuation coefficients to those of soft tissues

PHILIPS

B-mode Images of 3% Gelatin Phantoms with Different Oil Levels




- Higher attenuation for compressional wave in phantoms with higher oil concentration

PHILIPS

Shear Wave Displacement Profiles

One data acquisition at measurement depth 40mm





- Shear wave has higher frequency content and lower attenuation along the propagation path in 0% oil phantom than 20% oil phantom

PHILIPS

Shear Wave Phase Velocity as Function of Frequency and Voigt Model Fitting

One data acquisition at measurement depth 40mm
 3% gelatin + 0% oil 3% gelatin + 20% oil




- Shear wave speed increases more as a function of frequency (higher slope) in 20% oil phantom than 0% oil phantom


PHILIPS

Visco-elasticity Measurement for 3-7% Gelatin Concentration and 0 & 20% Oil

Shear modulus



Shear viscosity



- Shear modulus increases with gelatin percentage but decreases with castor oil percentage
- Shear viscosity increases with castor oil percentage

PHILIPS

Shear Wave Dispersion Measured by Shear Wave Speed Slope

Slope is measured by linear fitting shear wave phase velocity in the frequency range [70, 280]Hz

• 3% gelatin phantoms

Table 2. Shear modulus, shear viscosity and shear wave dispersion slope for 3% gelatin phantoms with four castor oil levels

	Castor oil level			
	0%	10%	20%	40%
Shear modulus (kPa)	2.01	1.68	1.10	0.88
Shear viscosity (Pa·s)	0.60	0.89	1.05	1.06
Dispersion slope (m/s/100 Hz)	0.06	0.15	0.27	0.32

• 6% gelatin phantoms

Table 3. Shear modulus, shear viscosity and shear wave dispersion slope for 6% gelatin phantoms with four castor oil levels

	Castor oil level			
	0%	10%	20%	40%
Shear modulus (kPa)	8.05	8.66	4.60	2.98
Shear viscosity (Pa·s)	1.14	2.25	2.50	1.97
Dispersion slope (m/s/100 Hz)	0.03	0.10	0.25	0.28

- Dispersion slope monotonically increases as a function of oil level
- 6% gelatin phantoms have dispersion slopes in a similar range as 3% gelatin phantom; but much higher absolute viscosity in 6% phantoms than in 3% phantoms
- Hypothesis: dispersion slope is a more independent measure of shear wave dispersion than shear viscosity derived from Voigt model fitting

PHILIPS

Conclusions and Discussions

- Tissue-mimicking viscous phantom were constructed using a gelatin base with additives including castor oil, alcohol (N-propanol), graphite powder and water
- Measured sound speeds and attenuation coefficients of fabricated phantoms were similar to those of soft tissues
- Reconstructed results show that gelatin concentration increases shear modulus, castor oil decreases shear modulus but increases shear wave dispersion and acoustic attenuation
- Study limitations
 - Phantom durability remains an issue (manufacturers such as CIRS may provide viscous phantoms of longer shelf life)
 - Shear wave frequency range characterized in this study is relatively narrow

PHILIPS

Main References on Phantom Development

- Hall TJ, Bilgen M, Insana MF, Krouskop TA. Phantom materials for elastography. IEEE Trans Ultrason Ferroelectr Freq Control 1997;44(6):1355-1365.
- Madsen EL, Zagzebski JA, Frank GR. Oil-in-gelatin dispersions for use as ultrasonically tissue-mimicking materials. Ultrasound Med Biol 1982;8(3):277-287.
- Barry CT, Mills B, Hah Z, Mooney RA, Ryan CK, Rubens DJ, Parker KJ. Shear wave dispersion measures liver steatosis. Ultrasound Med Biol 2012;38(2):175-182.
- Anderson PG, Rouze NC, Palmeri ML. Effect of graphite concentration on shear-wave speed in gelatin-based tissue-mimicking phantoms. Ultrason Imaging 2011;33(2):134-42.
- Zhang X, Qiang B, Greenleaf J. Comparison of the surface wave method and the indentation method for measuring the elasticity of gelatin phantoms of different concentrations. Ultrasonics 2011;51(2):157-164.

PHILIPS

